




MEDICO ELECTRODES

MEDICO ELECTRODES



Disposable ECG Electrodes

We Care For Your Heart !

HOLTER AND STRESS TEST ELECTRODES

A. Holter Electrodes



MLGLT-01, Liquid | MSGLT-01G, Solid
Size: 57 x 38 x 1 mm
Foam



MLGLT-03, Liquid | MSGLT-03G, Solid
Size: 50 x 48 x 1 mm
Foam



MLGLT-15, Liquid, | MSGLT-15G, Solid
Size: Ø 55 x 1 mm
Foam

B. Stress Electrodes



MLGLT-02, Liquid | MSGLT-04, Solid
Size: 50 x 38 x 1 mm | 50 x 35 x 1 mm
Foam



MSGLT-08G, Solid
Size: 45 x 42 x 1 mm
Foam




MSGLT-09, Solid
Size: 51 x 33 x 1 mm
Foam



MLGLT-14, Liquid, | MSGLT-14G, Solid
Size: Ø 50 x 1 mm
Foam

- Standard Packaging : 5 Electrodes/Card, 50 Electrodes/Pouch, 10 Pouches/Box, 12 Boxes/Carton.
- Customization Options Available.
- The Recommended usage indicated is only for guidance. It must be co-related with Patient's Physical Condition, Environmental Condition, Physician's Experience and Hospital Protocol.

All Product Sizes are mentioned in "hwxxt" format



For Illustration Purpose only

ALL PRODUCTS



Sensor
Ag+AgCl

RESTING & MONITORING ELECTRODES

A. Diagnostic Resting Electrodes



MUGST-64, Liquid | MSGST-64, Solid
Size: 44 x 32 x 1 mm
Foam, Banana Adaptor



MUGST-64M, Liquid | MSGST-64M, Solid
Size: 44 x 28 x 1 mm
Foam, Banana Adaptor



TAB ELECTRODES
(See Overleaf)

B. Monitoring Electrodes



MSGLT-05MG, Solid
Size: 40 x 36 x 1 mm
Foam



MSGST-06, Solid | MSGST-06M, Solid
Size: 44 x 32 x 1 mm | 44 x 28 x 1 mm
Foam



MSGLT-15DC, Solid
Size: Ø 55 mm
Non-Woven Cloth




MSGLT-17T, Solid
Size: 49 x 48 mm
Clear Tape

- Standard Packaging : 5 Electrodes/Card, 50 Electrodes/Pouch, 10 Pouches/Box, 12 Boxes/Carton.
- Customization Options Available.
- The Recommended usage indicated is only for guidance. It must be co-related with Patient's Physical Condition, Environmental Condition, Physician's Experience and Hospital Protocol.

All Product Sizes are mentioned in "hwxwt" format

SUBSTRATE OPTIONS



CONNECTOR OPTIONS



GEL OPTIONS



PEDIATRIC & NEO-NATAL ELECTRODES

A. Pediatric



MLGST-07, Liquid Foam
Size: 35 x 27 x 1 mm **P**



MSGST-37, Solid Foam
Size: 35 x 30 x 1 mm **P**




MSGST-27, Solid Foam
Size: 32 x 25 x 1 mm **P**

- P** Recommended for Pediatric Use
- N** Recommended for Neo-Natal Use


• Standard Packaging : 5 Electrodes/Card, 50 Electrodes/Pouch, 10 Pouches/Box, 12 Boxes/Cartron.

B. Neo-Natal - Prewired


Happy Trode



MPGBW - 30F, Hydrogel
Size: Ø 30 mm **N**




MPRYB - 78, Hydrogel Non-Woven Cloth
Size: 22 x 22 mm **N**




MPRYB - 79RTT, Hydrogel Tape Radiolucent
Size: 32 x 22 mm **N**



Happy Trode Plus⁺





MSRYB - 28RT, Hydrogel Non-Woven Cloth Radiolucent
Size: 28 x 20 mm **N**



MSGBW - 31RT, Hydrogel Non-Woven Cloth Radiolucent
Size: 19 x 16 mm **N**

- Standard Packaging : 3 Electrodes/Card, 3 Electrodes/Pouch, 100 Pouches/Box, 12 Boxes/Cartron.
- Customization Options Available.
- The Recommended usage indicated is only for guidance. It must be co-related with Patient's Physical Condition, Environmental Condition, Physician's Experience and Hospital Protocol.

All Product Sizes are mentioned in "hxwxh" format



TAB ELECTRODES

A. Wet Gel Tab Electrodes



TU-3728
Size: 37 x 28 mm
Foam
Radiolucent



- Standard Packaging : 10 Electrodes/Card, 100 Electrodes/Pouch, 20 Pouches/Box, 6 Boxes/Carton.

B. Solid Gel Tab Electrodes



3426MED
Size: 34 x 26 mm
Paper / PET
Radiolucent



3423MED
Size: 34 x 23 mm
Paper / PET
Radiolucent



3420MED | 3421MED | 3521MED
Size : 34 x 20 | 34 x 21 | 35 x 21 mm
Paper / PET
Radiolucent



3413MED
Size: 34 x 13 mm
Paper / PET
Radiolucent



- Standard Packaging : 10 Electrodes/Card, 100 Electrodes/Pouch, 10 Pouches/Box, 10 Boxes/Carton.
- Customization Options Available.
- The Recommended usage indicated is only for guidance. It must be co-related with Patient's Physical Condition, Environmental Condition, Physician's Experience and Hospital Protocol.

All Product Sizes are mentioned in "hxwxh" format

1 HYGIENE

Prevents Nosocomial Infection: Effectively eliminates the possible risks of cross contamination related to electrodes use



2 GOOD ADHESION

An effective adhesion throughout the ECG procedure



3 TRACE QUALITY

Skin impedance is effectively and immediately reduced by using Tab Electrodes, giving instantaneous and reliable signal pick up



4 PATIENT COMFORT

Safe and comfortable for all patient types



5 EASE OF USE

Easy to place and remove, optimising the ECG Procedure and improving clinical workflows




6 COST EFFECTIVE

Design is optimized to ensure cost effectiveness, saving valuable resources of the hospitals

OFF-SET ELECTRODES

A. Stress Electrodes



OTL-4227-S42 : Ø 42 mm
 OTL-4227-S50 : Ø 50 mm
 OTL-5235-S60 : Ø 60 mm
 Foam

P Recommended for Pediatric Use

• Standard Packaging : 1 Electrode/Card, 50 Electrodes/Pouch, 20 Pouches/Box, 6 Boxes/Carton

B. Resting Monitoring Electrodes



OTL-5235-S : Stainless Steel Stud
 OTL-5235 : Tab Connector
 OTL-5235-S4 : Banana Adaptor
 Size : 52 x 35 mm
 Foam



OTL-4227-S : Stainless Steel Stud
 OTL-4227 : Tab Connector
 OTL-4227-S4 : Banana Adaptor
 Size : 42 x 27 mm
 Foam

P

• Standard Packaging : 5 Electrodes/Card, 60 Electrodes/Pouch, 20 Pouches/Box, 6 Boxes/Carton

• Customization Options Available

• The Recommended usage indicated is only for guidance. It must be co-related with Patient's Physical Condition, Environmental Condition, Physician's Experience and Hospital Protocol.

All Product Sizes are mentioned in "h x w x t" format

1 Unique Off-Set connector reduces motion artifacts

Comfortable foam backing with enhanced adhesion ensures optimal fixation throughout the procedure **2**

3 Highly conductive wet gel ensures reliable traces

Also available in Radiolucent Connector **4**

ALL PRODUCTS




Sensor Ag+AgCl

SPECIALIZED ELECTRODES

A. Tab & Snap Electrodes



TSL-4532
Size: 45 x 32 mm
Foam
Radiolucent




- Standard Packaging : 5 Electrodes/Card, 60 Electrodes/Pouch, 20 Pouches/Box, 6 Boxes/Carton

B. Repositionable Electrodes



REP4032RT, Hydrogel
Size: 40 x 32 mm
Non-Woven Cloth
Radiolucent



- Standard Packaging : 3 Electrodes/Card, 3 Electrodes/Pouch, 200 Pouches/Box, 12 Boxes/Carton

All Product Sizes are mentioned in "h x w x t" format

- Customization Options Available
- The Recommended usage indicated is only for guidance. It must be co-related with Patient's Physical Condition, Environmental Condition, Physician's Experience and Hospital Protocol.



MEDICO ELECTRODES INTERNATIONAL LIMITED

Works: 142A/11&12, Noida Special Economic Zone, Noida 201305, India
Sales@MedicoElectrodes.com | www.MedicoElectrodes.com

