





MEDICO ELECTRODES

FlexStim[®]

Neurostimulation Electrodes



Say NO to ~~PAIN~~
Say YES to FLEXSTIM !



British-made Hydrogel



S5050C2SC

Size : 50 x 50 mm

Backing Material: Non-woven Cloth	4 pcs per Liner Card
Conductive Carbon film	4 pcs per Zip-lock Pouch
Gel: Conductive Hydrogel	500 pcs per Inner Box
Connector: Stainless Steel Stud	3000 pcs per Master Carton

S5050C2WC

Size : 50 x 50 mm

Backing Material: Non-woven Cloth	4 pcs per Liner Card
Conductive Carbon film	4 pcs per Zip-lock Pouch
Gel: Conductive Hydrogel	500 pcs per Inner Box
Connector: Lead Wire with 2 mm DIN	3000 pcs per Master Carton



C50C2SC

Size : Ø 50 mm

Backing Material: Non-woven Cloth	4 pcs per Liner Card
Conductive Carbon film	4 pcs per Zip-lock Pouch
Gel: Conductive Hydrogel	500 pcs per Inner Box
Connector: Stainless Steel Stud	3000 pcs per Master Carton

C50C2WC

Size : Ø 50 mm

Backing Material: Non-woven Cloth	4 pcs per Liner Card
Conductive Carbon film	4 pcs per Zip-lock Pouch
Gel: Conductive Hydrogel	500 pcs per Inner Box
Connector: Lead Wire with 2 mm DIN	3000 pcs per Master Carton



H5076C2SC

Size : 50 x 76 mm

Backing Material: Non-woven Cloth	1 pc per Liner Card
Conductive Carbon film	4 pcs per Zip-lock Pouch
Gel: Conductive Hydrogel	400 pcs per Inner Box
Connector: Stainless Steel Stud	2400 pcs per Master Carton

- Bag & Box Packaging Available
- Customization Options Available.
- The Recommended usage indicated is only for guidance. It must be co-related with Patient's Physical Condition, Environmental Condition, Physician's Experience and Hospital Protocol.

All Product Sizes are mentioned in "w x h" Format

1 Premium Quality: Designed for comfort and durability.

Durability: Reusable up to 20 times and more. **2**

3 British-Made Hydrogel: Premium Quality Hydrogel, manufactured at Medico's own plant in the United Kingdom.

Variety : Available in various shapes & sizes and with Silver coated carbon too. **4**

FDA 510(k)

ISO 13485

CE



Biocompatible

Re-usable

All Graphics are for Illustration purpose only



British-made Hydrogel



R50100C2SC

Size : 50 x 100 mm

Backing Material: Non-woven Cloth	2 pcs per Liner Card
Conductive Carbon film	4 pcs per Zip-lock Pouch
Gel: Conductive Hydrogel	400 pcs per Inner Box
Connector: Stainless Steel Stud	2400 pcs per Master Carton

R50100C2WC

Size : 50 x 100 mm

Backing Material: Non-woven Cloth	2 pcs per Liner Card
Conductive Carbon film	4 pcs per Zip-lock Pouch
Gel: Conductive Hydrogel	400 pcs per Inner Box
Connector: Lead Wire with 2 mm DIN	2400 pcs per Master Carton



O5090C2SC


Size : 50 x 90 mm

Backing Material: Non-woven Cloth	2 pcs per Liner Card
Conductive Carbon film	4 pcs per Zip-lock Pouch
Gel: Conductive Hydrogel	400 pcs per Inner Box
Connector: Stainless Steel Stud	2400 pcs per Master Carton

O5090C2WC

Size : 50 x 90 mm

Backing Material: Non-woven Cloth	2 pcs per Liner Card
Conductive Carbon film	4 pcs per Zip-lock Pouch
Gel: Conductive Hydrogel	400 pcs per Inner Box
Connector: Lead Wire with 2 mm DIN	2400 pcs per Master Carton




- Bag & Box Packaging Available
- Customization Options Available.
- The Recommended usage indicated is only for guidance. It must be co-related with Patient's Physical Condition, Environmental Condition, Physician's Experience and Hospital Protocol.

All Product Sizes are mentioned in "w x h" Format


SUBSTRATE OPTIONS



OEM Customization



Foam



Non-Woven Cloth

CONNECTOR OPTIONS



Carbon Pigtail wire
2mm DIN



Stainless Steel
Stud

CONDUCTOR OPTIONS



Carbon



Silver Coated Carbon

All Graphics are for Illustration purpose only



MEDICO ELECTRODES INTERNATIONAL LIMITED

Works: 142A/11&12, Noida Special Economic Zone, Noida 201305, India
Sales@MedicoElectrodes.com | www.MedicoElectrodes.com

