




# Medico Electrosurgical Neutral Electrodes

Grounding Plates (Universal Neutral Electrodes with Foam Backing)

Electrosurgical Electrodes having skin-friendly soft conductive adhesive hydrogel designed to improve contact quality thus ensuring a very low heat build-up on the skin!



# Medico Electrosurgical Neutral Electrodes

## Grounding Plates (Universal Neutral Electrodes with Foam Backing)

- ★ Medico Electrosurgical Electrodes are for single use and come with soft, flexible non-permeable polyethylene foam backing. The skin-friendly soft conductive adhesive hydrogel designed to improve contact quality ensures a very low heat build-up on the skin.
- ★ Medico Electrosurgical Electrodes are for universal applications i.e. for both adults and pediatric patients and come integrated with or without pre-attached lead wire.





### Key Benefits

Not made with natural rubber latex.

Universal in use (Both adult & pediatric patients).

Skin-friendly conductive hydrogel adhesive.

Shelf Life of 24 months from the date of production.

Product Picture	Product Item Code	Description	Packing
	<b>2820NT</b>	<p><b>Type:</b> Standard Electrosurgical Neutral Electrodes.  <b>Usage:</b> Universal (Both adult &amp; pediatric)  <b>Product Size:</b> 128 x 120 mm (height x width)  <b>Adhesive area:</b> 15,360 mm<sup>2</sup>  <b>Conductive area:</b> 10,200 mm<sup>2</sup></p>	<p>1 pc / Liner            10 pcs / Pouch            50 pcs / Inner carton            300 pcs / Master carton</p>
	<b>MNE-US</b>	<p><b>Type:</b> Split Electrosurgical Neutral Electrodes.  <b>Usage:</b> Universal (Both adult &amp; pediatric)  <b>Product Size:</b> 160 x 128 mm (height x width)  <b>Adhesive area:</b> 17,734 mm<sup>2</sup>  <b>Conductive area:</b> 10,807 mm<sup>2</sup></p>	<p>1 pc / Liner            10 pcs / Pouch            50 pcs / Inner carton            300 pcs / Master carton</p>
	<b>2820NW</b>	<p><b>Type:</b> Standard Electrosurgical Neutral Electrodes having integrated 3 meters wire with a standard Molex connector.  <b>Usage:</b> Universal (Both adult &amp; pediatric)  <b>Product Size:</b> 128 x 120 mm (height x width)  <b>Adhesive area:</b> 15,360 mm<sup>2</sup>  <b>Conductive area:</b> 10,200 mm<sup>2</sup></p>	<p>1 pc / Liner            1 pcs / Pouch            50 pcs / Inner carton            200 pcs / Master carton</p>
	<b>MNE-USW</b>	<p><b>Type:</b> Split Electrosurgical Neutral Electrodes having integrated 3 meters wire with a REM connector.  <b>Usage:</b> Universal (Both adult &amp; pediatric)  <b>Product Size:</b> 160 x 128 mm (height x width)  <b>Adhesive area:</b> 17,734 mm<sup>2</sup>  <b>Conductive area:</b> 10,807 mm<sup>2</sup></p>	<p>1 pc / Liner            1 pcs / Pouch            50 pcs / Inner carton            200 pcs / Master carton</p>



**MEDICO ELECTRODES INTERNATIONAL LTD.**

An ISO 13485 Company

Sales@MedicoElectrodes.com, www.MedicoElectrodes.com

