



FlexStim[®]

Neurostimulation Electrodes **TENS & EMS Electrodes**

- ★ Provides optimal & conductive interface between the TENS/EMS generator & Patient's skin.
 - ★ Strong adhesion & higher number of reusability.
 - ★ Available in skin friendly, breathable & flexible non-woven cloth for better skin conformability.
 - ★ Comes in various shapes & sizes.
 - ★ Foam Backing is also available.
-



MEDICO ELECTRODES



Product Specifications

Product Picture	Reference Code	Dimension	Components/Material	Standard Packaging
	C50C2SC	Ø 51 mm	Backing Material: Non-woven Cloth Conductive Carbon film Gel: Conductive Hydrogel Connector: Stainless Steel Snap	4 pcs per Liner Card. 4 pcs per Zip-lock Pouch. 60 pcs per Inner Box. 1800 pcs per Master Carton.
	C50C2WC	Ø 51 mm	Backing Material: Non-woven Cloth Conductive Carbon film Gel: Conductive Hydrogel Connector: Lead Wire with 2 mm DIN	
	S5050C2SC	51 x 51 mm	Backing Material: Non-woven Cloth Conductive Carbon film Gel: Conductive Hydrogel Connector: Stainless Steel Snap	
	S5050C2WC	51 x 51 mm	Backing Material: Non-woven Cloth Conductive Carbon film Gel: Conductive Hydrogel Connector: Lead Wire with 2 mm DIN	
	R50100C2SC	50 x 100 mm	Backing Material: Non-woven Cloth Conductive Carbon film Gel: Conductive Hydrogel Connector: Stainless Steel Snap	2 pcs per Liner Card. 4 pcs per Zip-lock Pouch. 40 pcs per Inner Box. 1200 pcs per Master Carton.
	R50100C2WC	50 x 100 mm	Backing Material: Non-woven Cloth Conductive Carbon film Gel: Conductive Hydrogel Connector: Lead Wire with 2 mm DIN	
	O5090C2SC	50 x 90 mm	Backing Material: Non-woven Cloth Conductive Carbon film Gel: Conductive Hydrogel Connector: Stainless Steel Snap	
	O5090C2WC	50 x 90 mm	Backing Material: Non-woven Cloth Conductive Carbon film Gel: Conductive Hydrogel Connector: Lead Wire with 2 mm DIN	
	H5076C2SC	50 x 76 mm	Backing Material: Non-woven Cloth Conductive Carbon film Gel: Conductive Hydrogel Connector: Stainless Steel Snap	1 pc per Liner Card. 2 pcs per Zip-lock Pouch. 50 pcs per Inner Box. 1500 pcs per Master Carton.

V100518

Available with Silver Coated Carbon Film also.

Packaging Customization options available.



MEDICO ELECTRODES INTERNATIONAL LTD.

An ISO 13485 Company

Sales@MedicoElectrodes.com www.MedicoElectrodes.com



MEDICO ELECTRODES