





Medico Electrosurgical Neutral Electrodes

Grounding Plates (Universal Neutral Electrodes with Foam Backing)

Electrosurgical Electrodes having skin-friendly soft conductive adhesive hydrogel designed to improve contact quality thus ensuring a very low heat build-up on the skin!



Medico Electrosurgical Neutral Electrodes

Grounding Plates (Universal Neutral Electrodes with Foam Backing)

- ★ Medico Electrosurgical Electrodes are for single use and come with soft, flexible non-permeable polyethylene foam backing. The skin-friendly soft conductive adhesive hydrogel designed to improve contact quality ensures a very low heat build-up on the skin.
- ★ Medico Electrosurgical Electrodes are for universal applications i.e. for both adults and pediatric patients and come integrated with or without pre-attached lead wire.





Key Benefits

Not made with natural rubber latex.

Universal in use (Both adult & pediatric patients).

Skin-friendly conductive hydrogel adhesive.

Shelf Life of 24 months from the date of production.

Product Picture	Product Item Code	Description	Packing
	2820NT	Type: Standard Electrosurgical Neutral Electrodes. Usage: Universal (Both adult & pediatric) Product Size: 128 x 120 mm (height x width) Adhesive area: 15,360 mm ² Conductive area: 10,200 mm ²	1 pc / Liner 10 pcs / Pouch 50 pcs / Inner carton 300 pcs / Master carton
	MNE-US	Type: Split Electrosurgical Neutral Electrodes. Usage: Universal (Both adult & pediatric) Product Size: 160 x 128 mm (height x width) Adhesive area: 17,734 mm ² Conductive area: 10,807 mm ²	1 pc / Liner 10 pcs / Pouch 50 pcs / Inner carton 300 pcs / Master carton
	2820NU	Type: Standard Electrosurgical Neutral Electrodes having integrated 3 meters wire with a standard Molex connector. Usage: Universal (Both adult & pediatric) Product Size: 128 x 120 mm (height x width) Adhesive area: 15,360 mm ² Conductive area: 10,200 mm ²	1 pc / Liner 1 pcs / Pouch 50 pcs / Inner carton 200 pcs / Master carton
	MNE-USW	Type: Split Electrosurgical Neutral Electrodes having integrated 3 meters wire with a REM connector. Usage: Universal (Both adult & pediatric) Product Size: 160 x 128 mm (height x width) Adhesive area: 17,734 mm ² Conductive area: 10,807 mm ²	1 pc / Liner 1 pcs / Pouch 50 pcs / Inner carton 200 pcs / Master carton

VZ51018