




# MEDICO ELECTRODES




## Disposable ECG Electrodes

We Care For Your Heart !

# HOLTER AND STRESS TEST ELECTRODES

## A. Holter Electrodes



**MLGLT-01, Liquid | MSGLT-01G, Solid**  
Size: 57 x 38 x 1 mm  
Foam



**MLGLT-03, Liquid | MSGLT-03G, Solid**  
Size: 50 x 48 x 1 mm  
Foam



**MLGLT-15, Liquid | MSGLT-15G, Solid**  
Size: Ø 55 x 1 mm  
Foam

## B. Stress Electrodes



**MLGLT-02, Liquid | MSGLT-04, Solid**  
Size: 50 x 38 x 1 mm | 50 x 35 x 1 mm  
Foam



**MSGLT-08G, Solid**  
Size: 45 x 42 x 1 mm  
Foam




**MSGLT-09, Solid**  
Size: 51 x 33 x 1 mm  
Foam



**MLGLT-14, Liquid | MSGLT-14G, Solid**  
Size: Ø 50 x 1 mm  
Foam

- Standard Packaging : 5 Electrodes/Card, 50 Electrodes/Pouch, 10 Pouches/Box, 12 Boxes/Carton.
- Customization Options Available.
- The Recommended usage indicated is only for guidance. It must be co-related with Patient's Physical Condition, Environmental Condition, Physician's Experience and Hospital Protocol.

All Product Sizes are mentioned in "h x w x t" format



For Illustration Purpose only

**ALL PRODUCTS**



Sensor  
Ag+AgCl

# RESTING & MONITORING ELECTRODES


## A. Diagnostic Resting Electrodes



**MLGST-64, Liquid | MSGST-64, Solid**  
Size: 44 x 32 x 1 mm  
Foam, Banana Adaptor



**MLGST-64M, Liquid | MSGST-64M, Solid**  
Size: 44 x 28 x 1 mm  
Foam, Banana Adaptor



**TAB ELECTRODES**  
(See Overleaf)

## B. Monitoring Electrodes



**MSGLT-05MG, Solid**  
Size: 40 x 36 x 1 mm  
Foam



**MSGST-06, Solid | MSGST-06M, Solid**  
Size: 44 x 32 x 1 mm | 44 x 28 x 1 mm  
Foam



**MSGLT-15DC, Solid**  
Size: Ø 55 mm  
Non-Woven Cloth




**MSGLT-17T, Solid**  
Size: 49 x 48 mm  
Clear Tape

- Standard Packaging : 5 Electrodes/Card, 50 Electrodes/Pouch, 10 Pouches/Box, 12 Boxes/Carton.
- Customization Options Available.
- The Recommended usage indicated is only for guidance. It must be co-related with Patient's Physical Condition, Environmental Condition, Physician's Experience and Hospital Protocol.

All Product Sizes are mentioned in "hxwxxt" format


### SUBSTRATE OPTIONS



### CONNECTOR OPTIONS



### GEL OPTIONS



# PEDIATRIC & NEO-NATAL ELECTRODES

## A. Pediatric



**MLGST-07, Liquid Foam**  
Size: 35 x 27 x 1 mm  
**P**



**MSGST-37, Solid Foam**  
Size: 35 x 30 x 1 mm  
**P**




**MSGST-27, Solid Foam**  
Size: 32 x 25 x 1 mm  
**P**

- P** Recommended for Pediatric Use
- N** Recommended for Neo-Natal Use


- Standard Packaging : 5 Electrodes/Card, 50 Electrodes/Pouch, 10 Pouches/Box, 12 Boxes/Carton.

## B. Neo-Natal - Prewired


### Happy Trode



**MPGBW - 30F, Hydrogel Foam**  
Size: Ø 30 mm  
**N**




**MPRYB - 78, Hydrogel Non-Woven Cloth**  
Size: 22 x 22 mm  
**N**




**MPRYB - 79RTT, Hydrogel Tape Radiolucent**  
Size: 32 x 22 mm  
**N**



### Happy Trode Plus<sup>+</sup>





**MSRYB - 28RT, Hydrogel Non-Woven Cloth Radiolucent**  
Size: 28 x 20 mm  
**N**



**MSGBW - 31RT, Hydrogel Non-Woven Cloth Radiolucent**  
Size: 19 x 16 mm  
**N**


- Standard Packaging : 3 Electrodes/Card, 3 Electrodes/Pouch, 100 Pouches/Box, 12 Boxes/Carton.
- Customization Options Available.
- The Recommended usage indicated is only for guidance. It must be co-related with Patient's Physical Condition, Environmental Condition, Physician's Experience and Hospital Protocol.

All Product Sizes are mentioned in "h x w x t" format




# TAB ELECTRODES

## A. Wet Gel Tab Electrodes




**TL-3728**  
**Size: 37 x 28 mm**  
**Foam**  
**Radiolucent**




- Standard Packaging : 10 Electrodes/Card, 100 Electrodes/Pouch, 20 Pouches/Box, 6 Boxes/Carton.


## B. Solid Gel Tab Electrodes



**3426MED**  
**Size: 34 x 26 mm**  
**Paper / PET**  
**Radiolucent**



**3423MED**  
**Size: 34 x 23 mm**  
**Paper / PET**  
**Radiolucent**



**3420MED | 3421MED | 3521MED**  
**Size : 34 x 20 | 34 x 21 | 35 x 21 mm**  
**Paper / PET**  
**Radiolucent**



**3413MED**  
**Size: 34 x 13 mm**  
**Paper / PET**  
**Radiolucent**



- Standard Packaging : 10 Electrodes/Card, 100 Electrodes/Pouch, 10 Pouches/Box, 10 Boxes/Carton.
- Customization Options Available.
- The Recommended usage indicated is only for guidance. It must be co-related with Patient's Physical Condition, Environmental Condition, Physician's Experience and Hospital Protocol.

All Product Sizes are mentioned in "hxwx" format

### 1 HYGIENE

Prevents Nosocomial Infection: Effectively eliminates the possible risks of cross contamination related to electrodes use



### 2 GOOD ADHESION

An effective adhesion throughout the ECG procedure



### 3 TRACE QUALITY

Skin impedance is effectively and immediately reduced by using Tab Electrodes, giving instantaneous and reliable signal pick up



### 4 PATIENT COMFORT

Safe and comfortable for all patient types



### 5 EASE OF USE

Easy to place and remove, optimising the ECG Procedure and improving clinical workflows




### 6 COST EFFECTIVE

Design is optimized to ensure cost effectiveness, saving valuable resources of the hospitals



# OFF-SET ELECTRODES

## A. Stress Electrodes



OTL-4227-S42 : Ø 42 mm  
 OTL-4227-S50 : Ø 50 mm  
 OTL-5235-S60 : Ø 60 mm  
 Foam


**P** Recommended for Pediatric Use

• Standard Packaging : 1 Electrode/Card, 50 Electrodes/Pouch, 20 Pouches/Box, 6 Boxes/Carton

## B. Resting Monitoring Electrodes



OTL-5235-S : Stainless Steel Stud  
 OTL-5235 : Tab Connector  
 OTL-5235-S4 : Banana Adaptor  
 Size : 52 x 35 mm  
 Foam



OTL-4227-S : Stainless Steel Stud  
 OTL-4227 : Tab Connector  
 OTL-4227-S4 : Banana Adaptor  
 Size : 42 x 27 mm  
 Foam


**P**

• Standard Packaging : 5 Electrodes/Card, 60 Electrodes/Pouch, 20 Pouches/Box, 6 Boxes/Carton

• Customization Options Available

• The Recommended usage indicated is only for guidance. It must be co-related with Patient's Physical Condition, Environmental Condition, Physician's Experience and Hospital Protocol.

All Product Sizes are mentioned in "hxwx" format



**ALL PRODUCTS**




Sensor Ag+AgCl

# SPECIALIZED ELECTRODES

## A. Tab & Snap Electrodes



**TSL-4532**  
 Size: 45 x 32 mm  
 Foam  
 Radiolucent




- Standard Packaging : 5 Electrodes/Card, 60 Electrodes/Pouch, 20 Pouches/Box, 6 Boxes/Carton

## B. Repositionable Electrodes



**REP4032RT, Hydrogel**  
 Size: 40 x 32 mm  
 Non-Woven Cloth  
 Radiolucent



- Standard Packaging : 3 Electrodes/Card, 3 Electrodes/Pouch, 200 Pouches/Box, 12 Boxes/Carton

All Product Sizes are mentioned in "hwxxt" format

- Customization Options Available
- The Recommended usage indicated is only for guidance. It must be co-related with Patient's Physical Condition, Environmental Condition, Physician's Experience and Hospital Protocol.



# MEDICO ELECTRODES INTERNATIONAL LIMITED

Works: 142A/11&12, Noida Special Economic Zone, Noida 201305, India  
Sales@MedicoElectrodes.com | [www.MedicoElectrodes.com](http://www.MedicoElectrodes.com)

