




FlexStim[®]

Neurostimulation Electrodes



Say NO to ~~PAIN~~
Say YES to FLEXSTIM !



British-made Hydrogel



S5050C2SC
Size : 50 x 50 mm

Backing Material: Non-woven Cloth | 4 pcs per Liner Card
Conductive Carbon film | 4 pcs per Zip-lock Pouch
Gel: Conductive Hydrogel | 500 pcs per Inner Box
Connector: Stainless Steel Stud | 3000 pcs per Master Carton



S5050C2WC
Size : 50 x 50 mm

Backing Material: Non-woven Cloth | 4 pcs per Liner Card
Conductive Carbon film | 4 pcs per Zip-lock Pouch
Gel: Conductive Hydrogel | 500 pcs per Inner Box
Connector: Lead Wire with 2 mm DIN | 3000 pcs per Master Carton



C50C2SC
Size : Ø 50 mm

Backing Material: Non-woven Cloth | 4 pcs per Liner Card
Conductive Carbon film | 4 pcs per Zip-lock Pouch
Gel: Conductive Hydrogel | 500 pcs per Inner Box
Connector: Stainless Steel Stud | 3000 pcs per Master Carton



C50C2WC
Size : Ø 50 mm

Backing Material: Non-woven Cloth | 4 pcs per Liner Card
Conductive Carbon film | 4 pcs per Zip-lock Pouch
Gel: Conductive Hydrogel | 500 pcs per Inner Box
Connector: Lead Wire with 2 mm DIN | 3000 pcs per Master Carton



H5076C2SC
Size : 50 x 76 mm

Backing Material: Non-woven Cloth | 1 pc per Liner Card
Conductive Carbon film | 4 pcs per Zip-lock Pouch
Gel: Conductive Hydrogel | 500 pcs per Inner Box
Connector: Stainless Steel Stud | 3000 pcs per Master Carton

- Bag & Box Packaging Available
- Customization Options Available.
- The Recommended usage indicated is only for guidance. It must be co-related with Patient's Physical Condition, Environmental Condition, Physician's Experience and Hospital Protocol.

All Product Sizes are mentioned in "w x h" Format



FDA 510(k)

ISO 13485



Biocompatible

Re-usable

All Graphics are for Illustration purpose only



British-made Hydrogel



R50100C2SC

Size : 50 x 100 mm

Backing Material: Non-woven Cloth	2 pcs per Liner Card
Conductive Carbon film	4 pcs per Zip-lock Pouch
Gel: Conductive Hydrogel	500 pcs per Inner Box
Connector: Stainless Steel Stud	3000 pcs per Master Carton

R50100C2WC

Size : 50 x 100 mm

Backing Material: Non-woven Cloth	2 pcs per Liner Card
Conductive Carbon film	4 pcs per Zip-lock Pouch
Gel: Conductive Hydrogel	500 pcs per Inner Box
Connector: Lead Wire with 2 mm DIN	3000 pcs per Master Carton



O5090C2SC


Size : 50 x 90 mm

Backing Material: Non-woven Cloth	2 pcs per Liner Card
Conductive Carbon film	4 pcs per Zip-lock Pouch
Gel: Conductive Hydrogel	500 pcs per Inner Box
Connector: Stainless Steel Stud	3000 pcs per Master Carton

O5090C2WC

Size : 50 x 90 mm


Backing Material: Non-woven Cloth	2 pcs per Liner Card
Conductive Carbon film	4 pcs per Zip-lock Pouch
Gel: Conductive Hydrogel	500 pcs per Inner Box
Connector: Lead Wire with 2 mm DIN	3000 pcs per Master Carton




- Bag & Box Packaging Available
- Customization Options Available.
- The Recommended usage indicated is only for guidance. It must be co-related with Patient's Physical Condition, Environmental Condition, Physician's Experience and Hospital Protocol.

All Product Sizes are mentioned in "w x h" Format


SUBSTRATE OPTIONS



OEM Customization



Foam



Non-Woven Cloth

CONNECTOR OPTIONS



Carbon Pigtail wire
2mm DIN



Stainless Steel
Stud

CONDUCTOR OPTIONS



Carbon



Silver Coated Carbon

All Graphics are for Illustration purpose only



MEDICO ELECTRODES INTERNATIONAL LIMITED

Works: 142A/11&12, Noida Special Economic Zone, Noida 201305, India
Sales@MedicoElectrodes.com | www.MedicoElectrodes.com

